

CODE 31497 MINI EAR SPECULUM diameter 4 mm - blue
APPLY THE FOLLOWING LABEL ON EVERY POLYBAG
SIZE OF LABEL: LENGTH105MM, HEIGHT52MM
IMPORTANT YOU WILL HAVE TO COMPLETE THE MANUFACTURING DATE: (XX-XXXX
MONTH-YEAR)



Conf 100 pz.
Bag of 100 pcs.



REF 31497

LOT 000037



MINI SPECULUM MONOUSO diam. 4 mm - blu
MINI EAR SPECULUM diameter 4 mm - blue
MINI-SPÉCULUM AURICULAIRE diam. 4 mm - bleu
MINI ESPÉCULO diám. 4 mm - azul




8 023279 314977


PE-LD PAP





Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com - Made in China





Rev. 3-11/2024