


Linux stethoscope



- Stetoscopio Linux
- Stéthoscope Linux
- Estetoscopio Linux
- Estetoscópio Linux
- Linux Stethoskop
- Linux Szimpla egyfejes fonendoszkóp
- Jednogłowicowy stetoskop Linux
- Linux Stetoskop med enkelt huvud
- Στηθοσκόπιο Linux
- Linux مسماع صدر



- Zinc alloy chestpiece with chrome plated silver finish
- Special design to get better acoustic effect pressing the back of head
- Padiglione in lega con finiture cromate
- Design studiato per ottimizzare l'acustica premendo sul padiglione

www.gimaitaly.com



Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

REF 32518 - 32524 - 32529

LOT Vedere etichetta bar code
See bar code label



Gima S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com
Made in China



S09111 - Rev.5 - 10/20
Size 16,3x32x43,5 cm



Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

32524 black/nero

32518 blue/blu

32529 red/rosso