

- 32728 Fluorimetro / Fluorimeter
 Adulti / Adults
 2 scale / 2 scales
 ○ 32734 Fluorimetro / Fluorimeter
 Adulti / Adults
 2 scale / 2 scales
 ○ 32709 Pediatrico / Pediatric
 Adulti / Adults
 2 scale / 2 scales

Tensionimètre Roma - Esfigmomanometro Roma
Roma Sphygmomanometer - Sfigmomanometro Roma

Roma Sphygmomanometer



Sfigmomanometro Roma
 Tensionimètre Roma
 Roma-Blutdruckmessgerät
 Esfigmomanometro Roma
 Esfigmomanometro Roma
 Плетизмограф Рима
 روميا ضغط الدم



- Lightweight, reliable, accurate
- Large dial Ø 60 mm
- Anti shock protective ring
- Screw and push button air release valve for one-hand use
- Wide range of models

- Leggere, affidabile, preciso
- Leggero, affidabile, preciso
- Facile da usare
- Facile da usare
- Valvola a pressione per uso a una mano
- Valvola a pressione per uso a una mano
- Leggere, affidabile, preciso
- Leggero, affidabile, preciso
- Facile da usare
- Facile da usare
- Valvola a pressione per uso a una mano
- Valvola a pressione per uso a una mano



www.gimaitaly.com



REF 32709 - 32727 - 32728
32734 - 32738

Gima S.p.A.
Via Marconi, 1 - 20080 Gessate (MI) Italy
Made in China



Via Marconi, 1 - 20080 Gessate (MI) Italy
gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

LOT Vedere etichetta bar code
See bar code label

2 Kids / Children
Table / Tableau

32734
Adulti / Adults
Adulti / Adults
Adulti / Adults



1 Table / Table

32727
Adulti / Adults
Adulti / Adults
Adulti / Adults



32709
Pediatrico / Pediatric
Pediatrico / Pediatric
Pediatrico / Pediatric

