

Sirio
Sphygmomanometer
SfignomanoMetro
Sirio
Tensiometre Sirio
Esfigmomanometro Sirio
مقياس ضغط الدم
Sirio-Blutdruckmessgerät
Πιεσομετρο Sirio

SfignomanoMetro Sirio



- Tensiometre Sirio
- Esfigmomanometro Sirio
- Πιεσομετρο Sirio



- Affidabile ed accurato
- Ampio quadrante per una facile lettura
- Bracciale in cotone con chiusura in velcro

www.gimaitaly.com




Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



Sirio sphygmomanometer




- Esfigmomanómetro Sirio
- Sirio-Blutdruckmessgerät
- مقياس ضغط الدم Sirio



- Reliable and accurate
- Large dial for easy reading
- Fabric cuff with velcro closure

www.gimaitaly.com



REF 32729

Assemblato da: / Assembled by:
Assemblé par: / Ensamblado por:
Gima S.p.A. - Gessate (MI) - Italy



Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

Componenti del kit / Kit components
Composants du kit / Componentes del kit
32904 MANOMETRO SIRIO CE 0476
32829 BRACCIALE COTONE + POLMONE 1 TUBO CE

