



GIMA

1 pz.
1 pc.



REF HS-70E (Gima 32744)

LOT 000005

SFIGMO DENVER SU CARRELLO

DENVER SPHYGMOMANOMETER - trolley - nylon cuff, TPU bladder

SPHYGMOMANOMÈTRE DENVER - sur socle à roulettes - brassard nylon, TPU vessie

ESFIGMO DENVER -sobre carro - brazal de nylon, cámara de air



HONSUN (NANTONG) Co., Ltd.
No.8 Tongxing Road Economic & Technological
Development Area, 226009 Nantong City, China
Made in China

CE 0123



Shanghai International Holding Corp.GmbH (Europe)
Eiffestrasse80,20537Hamburg, Germany

Rev.2.05.22

Importato da/Imported by/Importé de/Importado de
Gima S.p.A. - Via Marconi, 1
20060 - Gessate (MI) - Italy

