



www.gimaitaly.com

- Extra großer Quadrant Ø15 cm
- ABS-Korb zur Aufbewahrung der Manschette
- Kalibrierte Armmanschette aus Nylon
- Quadrante extra large Ø15 cm
- Estofo em ABS para bragaadeira
- Bragaadeira calibrada de nylon

PT

DE

Dallas-Blutdruckmessgerät - Wand Esfigmomanometro Dallas - muralha

Sfigmomanometro Dallas - da parede Tensiomètre Dallas - mural

IT

- Quadrante extra large Ø15 cm
- Cestello contenitore in ABS per bracciale
- Bracciale calibrato in nylon

FR

- Cadran extra large Ø15 cm
- Boîte en ABS pour le brassard
- Brassard étalonné en nylon



www.gimaitaly.com



Πιεσομετρο Dallas - τοίχου مقياس ضغط الدم Dallas - حائطي

GR

- Καντράν extra large Ø15 cm
- Καλάθι συσκευασίας με ABS για το περιβραχιόνιο
- Περιβραχιόνιο από νάιλον βαθμονομημένο

SA

- قرص مدرج واسع بقطر 15 سم
- وعاء حفظ للذراع من ABS
- ذراع مدرج من النايلون



www.gimaitaly.com



Dallas sphygmomanometer - wall type Esfigmomanómetro Dallas - mural

GB

- Extra large dial Ø15 cm
- ABS basket for cuff storage
- Nylon calibrated cuff

ES

- Pantalla muy grande Ø15 cm
- Contenedor en ABS para el brazalete
- Brazalete calibrado de nylon




www.gimaitaly.com



REF 32749

LOT Vedere etichetta
See label

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com
Made in China



Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com