

**GIMA**

conf. da 10 pz.
box of 10 pcs.
boîte de 10 pcs.

Rev.1.10.22



REF 5208 (GIMA 32858)

LOT 000006

PRESS SAN - GUAINA IGIENICA SFIGMOMANOMETRI
PRESS-SAN HYGIENIC PROTECTION OF SPHYGMOMANOMETER
PROTECTION HYGIÉNIQUE POUR TENSIOMÈTRE PRESS-SAN
PROTECCIÓN HIGIÉNICA para ESFIGMOMANOMETROS PRESS-SAN



Fabbricante/Manufacturer/Fabricant/Fabricante

Levetta s.a.s di Maffei Michele & C.

Via dei Cancellieri 6/8

51100 Pistoia (Italy)

Made in Italy

