

**REF** ECG300G (GIMA 33221)




ECG CONTEC 300G - 3 canali con display  
300G ECG - 3 channel with monitor  
ECG À 3 PISTES 300G avec moniteur  
ECG 300G - 3 canales con display



CONTEC MEDICAL SYSTEMS CO., LTD  
No. 112 Qinhuang West Street, Economic &  
Technical Development Zone, Qinhuangdao,  
Hebei Province, PEOPLE'S REPUBLIC OF CHINA

**EC** **REP**

Prolinx GmbH  
Brehmstr. 56, 40239, Duesseldorf,  
Germany  
Tel: 0049 21131054698  
E-mail: med@eulinx.eu



**CE** **0123**

Importato da/Imported by/Importé par/Importado por/  
GIMA S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) - Italy

**SN**

Vedere codice a barre  
See Bar Code  
Voir code à barres  
Ver código de barras  
Siehe Barcode



MADE IN CHINA

**SCATOLA**  
BOX



**RACCOLTA CARTA**  
PAPER COLLECTION

**SACCHETTO**  
**DI PLASTICA**  
PLASTIC BAG



**RACCOLTA PLASTICA**  
PLASTIC COLLECTION

**MANUALE UTENTE**  
USER MANUAL



**RACCOLTA CARTA**  
PAPER COLLECTION

Verifica le disposizioni del tuo comune  
per la gestione dei rifiuti  
Please refer to the provisions of your municipality  
concerning waste disposal