

REF BT005 (GIMA 33233)

LOT 000006



"P3" BATTERIA Li-ION RICARICABILE per 33232

"P3" RECHARGEABLE Li-ION BATTERY for 33232

"P3" BATTERIE RECHARGEABLE Li-ion pour 33232

"P3" BATERÍA Li-ION RECARGABLE para 33232



Importato da/Imported by/Importé par/Importado por:
GIMA S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) - Italy
Made in China

