

REF BIT0375 (Gima 33317)

LOT 000003

CAVO ECG per 33316 - ricambio
ECG CABLE for 33316 - spare
CÂBLE ECG pour 33316 - rechange
CABLE ECG para 33316 - recambio



Contec Medical Systems Co., Ltd
No.112 Qinhuang West Street, Economic & Technical
Development Zone, Qinhuangdao, Hebei Province,
066004, PEOPLE'S REPUBLIC OF CHINA
Made in China

EC **REP**

Shanghai International Holding Corp. GmbH (Europe)
Eiffestrasse 80, 20537 Hamburg, Germany



Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com



MD



55°C
-20°C

95%
0%

1060hPa
500hPa



Rev. 0-09/2023