

IT - ADULTI	Elettrodi con connettore per defibrillatori:
EN - ADULT	Electrodes with connector for defibrillators:
FR - ADULTE	Electrodes avec connecteur pour défibrillateurs:
DE - ERWACHSENE	Elektroden mit Konnektor für Defibrillatoren:
ES - ADULTO	Electrodos con conector para defibriladores:
PT - ADULTO	Electrodos com conector para desfibriladores:
NL - VOLWASSENEN	Plaat met connector voor defibrillators:
EL - ΕΝΗΛΙΚΑΣ	Ηλεκτρόδια με συνδετήρα για τους απινιδωτές:
RU - ДЛЯ ВЗРОСЛЫХ	Электроды для экг. дефибрилляции с коннектором
SV - VUXEN	Elektrodplattor med anslutning för defibrillator:


NIHON KOHDEN - WELCH ALLYN - INNOMED

DRAGER - S&W - MRL


I marchi, loghi e nomi commerciali menzionati sono proprietà dei rispettivi produttori
The mentioned brand, logos and commercial names are property of the corresponding manufacturers

REF F7955

LOT 1234XXXX



2024/08



(01) 1 8033003 26052 1 (17) 240800 (10) 1234XXXX



Via Costoli, 4 - 50039 - Vicchio
(Florence, Italy) - www.fiab.it

MADE IN ITALY


cod: 60010005
comp.: 60010005



Cod. 96502212ET9F-3° ed. - Compat.version: 48

NIHON KOHDEN - WELCH ALLYN - INNOMED

DRAGER - S&W - MRL



(01) 0 8033003 26052 4 (17) 240800 (10) 1234XXXX
cod.96503398ET9I
Compat. ver.: 48 F7955 comp.: 60010005 MADE IN ITALY