

mindray



Diagnostic Ultrasound System

Model: Z50

SN



XXX-XXXXXXXX



XXXX-XX-XX



047-030961-00(1.0)

100-240V ~ 50/60Hz 1.5-0.8A



SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD.

Mindray Building, Keji 12th Road South, High-Tech Industrial Park,
Nanshan, Shenzhen, 518057, P. R. China

EC REP

Shanghai International Holding Corp. GmbH(Europe)
Eiffestraße 80, 20537 Hamburg, Germany