

mindray



Digital Ultrasonic Diagnostic Imaging System

Model: XXXXXXXXX

SN



00000000



XXXXXXXXXX



047-027580-00(1.0)

100-240V ~ 50/60Hz 1.0-0.5A



SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD.

XXXXXXXXXX

EC REP

Shanghai International Holding Corp. GmbH(Europe)

Eiffestraße 80, 20537 Hamburg, Germany

XXXXXXXXXX