



1 pz.



REF 34332

LOT 000013

SET LARINGO "MAXLITE" F.O. PEDIATRICO 3 LAME MILLER 00,0,1 - 2,5V
MAXLITE F.O. PEDIATRIC SET 3 BLADES MILLER 00-0-1 - 2.5 V
SET F.O. MAXLITE PÉDIATRIQUE 3 lames: Miller 00-0-1 - 2,5V
SET LARINGO "MAXLITE" F.O. PEDIÁTRICO 3 HOJAS MILLER 00,0,1 - 2,5V



GIMA S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
Made in Pakistan - www.gimaitaly.com

