

Laringoscopi a Fibre Ottiche




- Laryngoscope à fibre optique
- Laringoscópios em fibra ótica
- Laryngoskopy światłowodowe
- Λαρυγγοσκόπιο οπτικής ίνας



- Luce fredda ad elevata luminosità - fibra Ø 4 mm
- Compatibile con tutti i sistemi F.O. "Green"
- Fibra facilmente rimovibile ed autoclavabile

www.gimaitaly.com




Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com
Made in Pakistan

VEDERE ETICHETTA BAR CODE
SEE BAR CODE LABEL



Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

Fiber Optic Laryngoscopes



- Laringoscopios de fibra óptica
- LWL-Laryngoskope
- Laringoscop cu fibră optică
- منظار الحنجرة بألياف بصرية



- Superior cold, bright light Ø 4 mm fibre
- Compatible with all "Green" F.O. systems
- Autoclavable, replaceable light guide

www.gimaitaly.com



Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

REF 34451 - 34452 - 34453
34454 - 34455

Fiber Optic Laryngoscopes - Laringoscopi a Fibre Ottiche

- 34451 Pediatric set 2 blades Miller 0-1
- 34452 Adult set 3 blades M.C. Intosh 2-3-4
- 34453 Adult set 4 blades M.C. Intosh 1-2-3-4
- 34454 Pediatric set 3 blades Miller 0-1 M.C. Intosh 0
- Doctor set 3 blades Miller 1-2-3

Size 19x23,5x4,8 cm
509037 - Rev. 5 - 04/20