



1 pz.



REF 34456

LOT 000008

MANICO IN METALLO RICARICABILE 3,5V - pediatrico
"GIMA GREEN" RE-CHARGEABLE HANDLE 3.5V - pediatric
MANCHE RECHARGEABLE "GIMA GREEN" pédiatrique (3,5V)
MANGO "GIMA GREEN" RECARGABLE 3,5V - pediátrico



GIMA S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) - Italy
Made in Pakistan - www.gimaitaly.com

