


REF 808-069-25 (GIMA 34466) **LOT** 000012



BATTERIA RICARICABILE PEDIATRICA Li-ion 2,5V
PEDIATRIC RE-CHARGEABLE Li-ion BATTERY 2.5V 1000mAh
BATTERIE RECHARGEABLE PÉDIATRIQUE Li-ion 2,5V 1000mAh
"P3" BATERÍA RECARGABLE PEDIÁTRICA Li-ion 2,5V 1000mAh

Fabbricante/Manufacturer/Fabricant/Fabricante:
MEDICAL DEVICES (PVT) LTD. Wazirabad Road, Ugoki, Sialkot Pakistan
info@mdevices.com

2.5 V 1000mAh
Tempo di ricarica: 2,5 ore
Charging time: 2.5 hours
Temps de charge : 2,5 heures
Tiempo de carga: 2,5 horas

Importato da/Imported by/ Importé par/ Importado por:
GIMA S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) - Italy – Made in Pakistan
gima@gimaitaly.com - export@gimaitaly.com - www.gimaitaly.com

