

REF 15040054 (GIMA 35097)

MODELLO/MODEL/MODÈLE/MODELO: KS-R02

SONDA SpO2 PEDIATRICA in GOMMA - riusabile per 35095
SpO2 PEDIATRIC PROBE - soft rubber, reusable - for 35095
CAPTEUR PÉDIATRIQUE EN CAOUTCHOUC - réutilisable pour 35095
SONDA SpO2 PEDIÁTRICA de GOMA SUAVERE - reutilizable para 35095



Shenzhen Creative Industry Co., Ltd.
Floor 5, BLD 9, BaiWangxin High-Tech Industrial Park,
Songbai Road, Xili Street, Nanshan District,
518110 Shenzhen, P.R. China
MADE IN CHINA

EC REP

Shanghai International Holding Corp. GmbH (Europe)
Eiffelstrasse 80, 20537 Hamburg, Germany

LOT 000023



Importato da/Imported by/Importé par/Importado por:
GIMA S.p.A. - Via Marconi, 1 - 20060 Gessate (MI) - Italy



BUSTINA
POUCH



RACCOLTA PLASTICA
PLASTIC COLLECTION

Verifica le disposizioni del tuo comune per la gestione dei rifiuti
Please refer to the provisions of your municipality
concerning waste disposal



GIMA
www.gimaitaly.com



C E 0123



IPX2