

REF 15040033 (GIMA 35099)

LOT 000022

GIMA
www.gimaitaly.com

MODELLO/MODEL/MODÈLE/MODELO: KS-YW02

SONDA SpO2 NEONATALE - riusabile per 35095

SpO2 NEONATAL PROBE - reusable - for 35095

CAPTEUR SpO2 NÉONATAL réutilisable pour 35095

SONDA SpO2 NEONATAL - reutilizable para 35095



Shenzhen Creative Industry Co., Ltd.

Floor 5, BLD 9, BaiWangxin High-Tech Industrial Park,

Songbai Road, Xili Street, Nanshan District,

518110 Shenzhen, P.R. China

MADE IN CHINA

EC **REP**

Shanghai International Holding Corp. GmbH (Europe)

Eiffestrasse 80, 20537 Hamburg, Germany



Importato da / Imported by / Importé par / Importado por:

Gima S.p.A. Via Marconi, 1 - 20060 Gessate (MI) Italy



CE 0123



IPX2



BUSTINA
POUCH



RACCOLTA PLASTICA
PLASTIC COLLECTION

ISTRUZIONI PER L'USO
USER INSTRUCTIONS



RACCOLTA CARTA
PAPER COLLECTION

Verifica le disposizioni del tuo comune per la gestione dei rifiuti
Please refer to the provisions of your municipality
concerning waste disposal



Rev.5-11/2023