

REF

21152

Quantità/Quantity

800



IVD

STERILE R

made
in Italy

Descrizione

**PROMED® PR. 16x100 mm 10 ml CILINDRICA IN POLISTIROLO CON TAPPO
ROSSO + ET. STERILE SFUSA**

Description

**PROMED® TEST TUBE 16x100 mm 10 ml CYLINDRICAL IN POLYSTYRENE WITH
RED CAP AND LABEL STERILE IN BULK**

LOT

XXXXXXXX



02/2029

Operatore/Operator



F.L. MEDICAL S.r.l. Unipersonale - Via Enrico Mattei, 20 - 35038 Torreglia (PD) - ITALY
Phone +390495211475 - Fax +390495212566 - mail info@flmedical.com