

REF

Quantità/Quantity



27142

45

made
in Italy

Descrizione

**PROMED® PORTAPROVETTE NON MONTATI IN POLIPROPILENE
PER PROVETTE Ø 12-13 mm A 90 POSTI, AUTOCLAVABILE, COLORE
VERDE**

Description

**PROMED® NOT-ASSEMBLED TEST TUBE RACK IN POLYPROPYLENE
FOR TUBE Ø 12-13 mm 90 PLACES, AUTOCLAVABLE, GREEN
COLOUR**

LOT

XXXXXXXX



Operatore/Operator



F.L. MEDICAL S.r.l. Unipersonale - Via Enrico Mattei, 20 - 35038 Torreglia (PD) - ITALY
Phone +390495211475 - Fax +390495212566 - mail info@flmedical.com

0