



GIMA

1 pz.
1 pc.



GIMA 44740

LOT 000005




LETTO A 3 SEZIONI AD ALTEZZA VARIABILE - 2 manovelle - con ruote
3 JOINTS VARIABLE HEIGHT PATIENT BED - 2 cranks - castors
LIT 3 ARTICULATIONS HAUTEUR VARIABLE - 2 crémaillères - sur roulettes
CAMA 3 ARTICULACIONES ALTURA VARIABLE - 2 manivelas - sobre



MD



Jiangsu Saikang Medical Equipment Co.,Ltd
No. 35 Lehong Road, Modern Agriculture
Demonstration Park, Zhangjiagang City,
Jiangsu Province, China
Made in China



Rev.2.01.25

REF V2K5C



Sungo Europe B.V.
Fascinatio Boulevard 522, Unit 1.7, 2909VA
Capelle aan den IJssel, The Netherlands



Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (Mi) - Italy
gima@gimaitaly.com -
export@gimaitaly.com
www.gimaitaly.com

