



GIMA

1 pz.

1 pc.

REF AP338 (GIMA 24486)

LOT Vedere etichetta / See label

"P16" ETILOMETRO PROFESSIONALE AP338 con stampante

"P16" AP338 PROFESSIONAL ALCOHOL TESTER EASY RECALIBRATION with printer

ÉTHYLOTEST PROFESSIONNEL AP338 CALIBRAGE FACILE avec imprimante

PROBADOR DE ALCOHOL PROFESIONAL AP338 RECALIBRACIÓN FÁCIL con impresora



Fabbricante/Manufacturer/Fabricant/Fabricante:

A&A Product Manufacturing Limited

Workshop 9, 18/F, Tower A, New Trade Plaza, 6 On Ping Street,
ShaTin, New Territories, Hong Kong

sales@aap.com.hk

Made in China

Importato da/Imported by/Importé par/Importado por:

Gima S.p.A. Via Marconi, 1 - 20060 Gessate (MI) Italy

gima@gimaitaly.com - export@gimaitaly.com

www.gimaitaly.com

