



**GIMA**

REF 25584  
LOT

Model:  
HTD8816C

FC



RoHS



SN



**Non-Contact Infrared Body Thermometer**

Accuracy: 35,0°C-42,0°C: ± 0,2°C / 95,0°F-107,6°F: ± 0,4°F

Power: DC3V (AAAx2)

Measuring site: Forehead

Range: 34,0-43,0°C (93,2°F-109,4°F)



CE

IP22

0598

**HeTaiDa Technology Co., Ltd**  
4F BaShiDa High-Tech Park, XiangDong Industrial Area,  
DaLingShan Town, DongGuan City, GuangDong, China  
Made in China

**WellKang Ltd.**  
16 Castle St, Dover, CT16 1PW Dover,  
CT16 1PW, UK