



GIMA

Bilancia meccanica/Mechanical scale

Capacità: 125 Kg Precisione: 1 Kg - Capacity: 125 Kg Division: 1Kg

REF 7720 (27230)

LOT



Manufacturer/Produttore/Fabricant/Fabricante:

Momert Co.Ltd. 2400 Dunaujváros - HUNGARY, Papírgyári út 12-14

Made in Hungary

Distribuito da/Distributed by/Distribué par/Distribuido por

Gima S.p.A. Via marconi, 1

20060 - Gessate (MI) - Italy

www.gimaitaly.com

Rev. 1.03.21