



GIMA



REF TCS-200-RT
(GIMA 27289)

LOT 000007



Capacity - Portata - Portée - Capacidad - Capacidade : 200 kg -
Sensitivity - Sensibilità - Sensibilité - Sensibilidad - Sensibilidade: 100 g

Fabbricante/Manufacturer/Fabricant/Fabricante/Fabricante

Wuxi Weigher Factory Co., Ltd.

Majiaqiao, Guangfeng, Guangyi – 214011 Wuxi, Jiangsu

PEOPLE'S REPUBLIC OF CHINA - scale8@aliyun.com

Made in China

Importato da/Imported by/Importé par/Importado por/Importado por

Gima S.p.A.

Via Marconi,1 - 20060 Gessate (Mi) – Italy

gima@gimaitaly.com - export@gimaitaly.com

www.gimaitaly.com

Rev.5.10.24