



GIMA

www.gimaitaly.com

GIMA 27708

Rev.4.03.25

LOT 000032



Max load: 75 Kg

Carico massimo: 75 Kg

Charge maximale: 75 Kg

Carga máxima: 75 Kg

MD



Foshan Dongfang Medical
Equipment Manufactory (LTD.)
5 Xingfu Road, Zone A, Shishan
Industrial Park, Nanhai District,
Foshan City, P.R. China, 528000
Made in China

REF FS980LA-35



EC

REP

MedNet ECREPGmbH

Borkstrasse 10

48163 Muenster. Germany

Importato da/Imported by

Importé par/Importado por

Gima S.p.A

Via Marconi, 1 – 20060 Gessate (Mi) – Italy

gima@gimaitaly.com - export@gimaitaly.com

www.gimaitaly.com