



GIMA

www.gimaitaly.com

GIMA 27709



LOT 000093

Max load:100Kg

Carico massimo:100Kg

Charge maximale:100 Kg

Carga máxima:100Kg



MD



Foshan Dongfang Medical
Equipment Manufactory (LTD.)
5 Xingfu Road, Zone A, Shishan
Industrial Park, Nanhai District,
Foshan City, P.R. China, 528000
Made in China



REF FS809-46



MedNet EC-Rep GmbH Borkstrasse 10
48163 Muenster. Germany

Importato da/Imported by/Importé par/Importado por
Gima S.p.A

Via Marconi, 1 – 20060 Gessate (Mi) – Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

Rev.6-04.25