



GIMA

www.gimaitaly.com

GIMA 27715

MD



LOT 000059



Max load: 100 Kg

Carico massimo: 100 Kg

Charge maximale: 100 Kg

Carga máxima: 100 Kg



Rev.5.06.25

Foshan Dongfang

Medical Equipment Manufactory (LTD.)

5 Xingfu Road, Zone A, S

hishan Industrial Park, Nanhai District,

Foshan City, P.R. China, 528000

Made in China

REF FS908AQ-46



EC

REP

MedNet EC-Rep GmbH

Borkstrasse 10 - 48163 Muenster. Germany

Importato da/Imported by/ Importé par/

Importado por

Gima S.p.A

Via Marconi, 1 – 20060 Gessate (Mi) Italy

gima@gimaitaly.com -

export@gimaitaly.com

www.gimaitaly.com