



GIMA

www.gimaitaly.com

GIMA 27716



LOT 000039

Max load: 125 Kg

Carico massimo: 125 Kg

Charge maximale: 125 Kg

Carga máxima: 125 Kg



MD



Rev.5.04.25

Foshan Dongfang

Medical Equipment Manufactory (LTD.)

5 Xingfu Road, Zone A, S

hishan Industrial Park, Nanhai District,

Foshan City, P.R. China, 528000

Made in China

REF FS951B-56



MedNet EC-Rep GmbH

Borkstrasse 10 - 48163 Muenster. Germany

Importato da/Imported by/

Importé par/Importado por

Gima S.p.A

Via Marconi, 1 – 20060 Gessate (Mi) Italy

gima@gimaitaly.com - export@gimaitaly.com

www.gimaitaly.com