



GIMA

www.gimaitaly.com

REF **FS957LQ**
(27720)

LOT **000013**



Rev.2.09.21



Max load: 100 Kg

Carico massimo: 100 Kg

Charge maximale: 100 Kg

Carga máxima: 100 Kg

MD



Foshan Dongfang Medical
Equipment Manufactory (LTD.)
5 Xingfu Road, Zone A, Shishan
Industrial Park, Nanhai District,
Foshan City, P.R. China, 528000
Made in China



MedNet GmbH Borkstrasse 10
48163 Muenster. Germany

Importato da/Imported by
Importé par/Importado por
Gima S.p.A



Via Marconi, 1 – 20060 Gessate (Mi) – Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com