



GIMA

www.gimaitaly.com

REF FS935S(L)
(GIMA 27799)

LOT 000010

MD



Max load: 100 Kg
Carico massimo: 100 Kg
Charge maximale: 100 Kg
Carga máxima: 100 Kg

Rev.3.02.22



Foshan Dongfang Medical
Equipment Manufactory (LTD.)
5 Xingfu Road, Zone A, Shishan
Industrial Park, Nanhai District,
Foshan City, P.R. China, 528000
Made in China



MedNet EC-Rep GmbH
Borkstrasse 10
48163 Muenster. Germany

Importato da/Imported by
Importé par/Importado por
Gima S.p.A

Via Marconi, 1 – 20060 Gessate (Mi) – Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com