



GIMA

REF


QDC-303-P4000IIE (B) (Gima 28566)

LOT


Rev.2.02.23

MATERASSO A BOLLE + COMPRESSORE (28564+28565)
BUBBLE AIR MATTRESS + COMPRESSOR (28564+28565)
MATELAS À BULLES + COMPRESSEUR (28564+28565)
COLCHONETA DE BURBUJAS + COMPRESOR (28564+28565)

SN



Guangdong Yuehua Medical Instrument Factory Co., Ltd.
Rongsheng Science and Technology Zone, Daxue Road, Shantou, Guangdong, China
Made in China



SUNGO Europe B.V. Fascinatio Boulevard 522, Unit 1.7, 2909VA Capelle aan den IJssel, The Netherlands

MD

Importato da/Imported by/Importé par/Importado por
Gima S.p.A Via Marconi, 1 – 20060 Gessate (Mi) – Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

POWER 220V AC 50Hz 0.1A

