



GIMA

REF

QDC-8010 + P3000A2QB3
(GIMA 28580)

LOT

000003

KIT 8010 MATERASSO A PRESSIONE ALTERNATA NYLON PVC+POMPA
KIT 8010 NYLON PVC ALTERNATING PRESSURE MATTRESS + PUMP
KIT 8010 MATELAS À PRESSION ALTERNÉE NYLON PVC + POMPE
KIT 8010 COLCHÓN PRESIÓN ALTERNADA NYLON PVC + BOMBA

POWER 220V AC 50Hz 0.1A



Guangdong Yuehua Medical Instrument Factory Co., Ltd.
Rongsheng Science and Technology Zone, Daxue Road,
Shantou, Guangdong, China
Made in China

Rev.2.03.25



Eunitor GmbH

Kennedydamm, 5, Düsseldorf, 40476, Germany

Importato da/Imported by/Importé par/Importado por
Gima S.p.A
Via Marconi, 1 – 20060 Gessate (Mi) – Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

