



GIMA

REF

QDC-8080+P3000A2QB3
(GIMA 28582)

LOT

000003

KIT 8080 MATERASSO A PRESSIONE ALTERNATA NYLON TPU+POMPA
KIT 8080 NYLON TPU ALTERNATING PRESSURE MATTRESS + PUMP
KIT 8080 MATELAS À PRESSION ALTERNÉE NYLON TPU + POMPE
KIT 8080 COLCHÓN PRESIÓN ALTERNADA NYLON TPU + BOMBA

POWER 220V AC 50Hz 0.1A



Guangdong Yuehua Medical Instrument Factory Co., Ltd.
Rongsheng Science and Technology Zone, Daxue Road,
Shantou, Guangdong, China
Made in China

Rev.1.02.23

MD

Eunitor GmbH

Kennedydamm, 5, Düsseldorf, 40476, Germany

Importato da/Imported by/Importé par/Importado por
Gima S.p.A
Via Marconi, 1 – 20060 Gessate (Mi) – Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

