



GIMA

Rev.4.01.23



REF 29777

LOT 000013

VETRINI MOLATI-SABBIATI su 2 lati 26x76x1,2 mm

SLIDES 26x76 mm - twin frosted in both ends

LAMES POUR MICROSCOPE 26x76 mm - Extrémités dépolies, bande sablée bords rodés

PORTAOBJETOS BISELADOS-PULIDOS los 2 lados 26x76 mm



Importato da/Imported by/Importé par/Importado
por Importado por/Εισαγωγή από/

Gima S.p.A.

Via Marconi, 1 - 20060 Gessate (MI) Italy

gima@gimaitaly.com - export@gimaitaly.com

www.gimaitaly.com

Made in China



8 023279 297775