


1.5mm



CONTEC

Pulse Oximeter

Model: CMS60D

Safety Class: Type BF 

Battery Voltage: DC 2.6 V-3.6 V Operating Current: ≤ 100 mA



Contec Medical Systems Co., Ltd.

IP22

No.112 Qinhuang West Street, Economic & Technical
Development Zone, Qinhuangdao, Hebei Province,
PEOPLE'S REPUBLIC OF CHINA

MD

EC REP

Prolinx GmbH

Brehmstr. 56, 40239 Duesseldorf Germany



SN

CE 0123

