



GIMA

www.gimaitaly.com

GIMA 42903

MD



LOT 000010

Max load: 100 Kg
Carico massimo: 100 Kg
Charge maximale: 100 Kg
Carga máxima: 100 Kg



Rev.4.06.25

Foshan Dongfang
Medical Equipment Manufactory (LTD.)
5 Xingfu Road, Zone A, S
hishan Industrial Park, Nanhai District,
Foshan City, P.R. China, 528000
Made in China

REF FS7972L



MedNet EC-Rep GmbH
Borkstrasse 10 - 48163 Muenster. Germany
Importato da/Imported by/ Importé
par/Importado por
Gima S.p.A
Via Marconi, 1 – 20060 Gessate (Mi) Italy
gima@gimaitaly.com -
export@gimaitaly.com
www.gimaitaly.com