



GIMA

www.gimaitaly.com

REF FS7981L (GIMA 42906)

LOT 000010



Max load: 150 Kg
Carico massimo: 150 Kg
Charge maximale: 150 Kg
Carga máxima: 150 Kg



Rev.2.12.22



MD

Foshan Dongfang
Medical Equipment Manufactory (LTD.)
5 Xingfu Road, Zone A, S
hishan Industrial Park, Nanhai District,
Foshan City, P.R. China, 528000
Made in China



MedNet EC-Rep GmbH
Borkstrasse 10 - 48163 Muenster. Germany

Importato da/Imported by/
Importé par/Importado por
Gima S.p.A

Via Marconi, 1 – 20060 Gessate (Mi)
Italy

gima@gimaitaly.com -
export@gimaitaly.com
www.gimaitaly.com