



GIMA

www.gimaitaly.com

GIMA 43070

LOT 000016

MD



Max load: 100 Kg

Carico massimo: 100 Kg

Charge maximale: 100 Kg

Carga máxima: 100 Kg



Foshan Dongfang Medical
Equipment Manufactory (LTD.)
5 Xingfu Road, Zone A, Shishan
Industrial Park, Nanhai District,
Foshan City, P.R. China, 528000
Made in China

REF FS922



EC

REP

MedNet EC-Rep GmbH
Borkstrasse 10 - 48163
Muenster. Germany

Rev.5.06.25

Importato da/Imported by
Importé par/Importado por
Gima S.p.A

Via Marconi, 1 – 20060 Gessate (Mi) – Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com