



GIMA

www.gimaitaly.com

GIMA 43225

LOT 000011

MD

CE



Max load: 100 Kg

Carico massimo: 100 Kg

Charge maximale: 100 Kg

Carga máxima: 100 Kg



Rev.3.03.25



Foshan Dongfang Medical
Equipment Manufactory (LTD.)
5 Xingfu Road, Zone A, Shishan
Industrial Park, Nanhai District,
Foshan City, P.R. China, 528000
Made in China

REF FS863ABJ-46

EC

REP

MedNet EC-Rep GmbH
Borkstrasse 10
48163 Muenster. Germany

Importato da/Imported by
Importé par/Importado por
Gima S.p.A

Via Marconi, 1 – 20060 Gessate (Mi) – Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com