



GIMA



www.gimaitaly.com

Rev 0.12.22

REF YK9063 (GIMA 43235)



MD



LOT 000003

Max load/Carico massimo/ Charge maximale/Carga máxima: 120 Kg



INTCO MEDICAL (HK) CO., LTD
FLAT/RM 19C, LOCKHART CENTRE, 301-307 LOCKHART ROAD,
WAN CHAI, HONG KONG, P.R. China
Made in China



MedNet EC-REP GmbH
Borkstrasse 10, 48163 Muenster, Germany

Importato da/Imported by/Importé par/Importado por:

Gima S.p.A.

Via Marconi, 1 - 20060 Gessate (MI) Italy

gima@gimaitaly.com - export@gimaitaly.com

www.gimaitaly.com